

## TURBINE HANDPIECES

synea VISION



The new range of Synea Vision turbine handpieces: Innovative technology and a comprehensive range of products for the highest demands.

- > 100% shadowless illumination of the treatment site with the ring LED (with the RotoQuick system only)
- > Optimum cooling of the treatment site with 5-fold spray
- > Four head sizes for every application case
- > Scratch-resistant surface improves aesthetics even more – and extends their life
- > Fatigue-free work
- > Optimum service life and unexpected quietness during operation with ceramic ball bearings
- > Thermo washer disinfectable and sterilizable
- > 24 months warranty

### SYNEA VISION TURBINE HANDPIECES



Ring LED+,  
5-fold spray, 22 W,  
head size Ø 13mm,  
press-button chuck system for  
FG burs Ø 1,6 mm

**TK-100 L**  
REF 30021000



Ring LED+,  
5-fold spray, 21 W,  
head size Ø 11,5mm,  
press-button chuck system for  
FG burs Ø 1,6 mm

**TK-98 L**  
REF 30019000



Ring LED+,  
5-fold spray, 18 W,  
head size Ø 10 mm,  
press-button chuck system for  
FG burs Ø 1,6 mm

**TK-97 L**  
REF 30017000



LED+,  
5-fold spray, 12 W,  
head size Ø 9 mm,  
press-button chuck system for  
FG burs Ø 1,6 mm

**TK-94 L**  
REF 30016000



See Roto Quick couplings

### SYNEA VISION TURBINE HANDPIECES WITH MULTIFLEX® CONNECTION



glass rod, 5-fold spray,  
22 W, head size Ø 13 mm,  
press-button chuck system for  
FG turbine burs Ø 1,6 mm

**TK-100 LM**  
REF 30022000



glass rod, 5-fold spray,  
21 W, head size Ø 11,5 mm,  
press-button chuck system for  
FG turbine burs Ø 1,6 mm

**TK-98 LM**  
REF 30020000



glass rod, 5-fold spray,  
18 W, head size Ø 10 mm,  
press-button chuck system for  
FG turbine burs Ø 1,6 mm

**TK-97 LM**  
REF 30018000



glass rod, 5-fold spray,  
12 W, head size Ø 9 mm,  
press-button chuck system for  
FG turbine burs Ø 1,6 mm

**TK-94 LM**  
REF 30031000



See W&H coupling with LED light

led+ with ring LED+ light

led+ with LED+ light



with light



without light



integrated spray



thermo washer  
disinfectable



sterilizable

## CONTRA-ANGLE AND STRAIGHT HANDPIECES


synea VISION



Synea Vision contra-angle and straight handpieces: reliable, ergonomic and outstanding quality.

- > Two head sizes for individual needs
- > Optimum service life and quiet operation due to ceramic ball bearings and innovative gearing technology (WK-99LT / WK-93/LT)
- > Easy to clean, comprehensive protection and long-term value retention thanks to new surface coating
- > Excellent grip and ergonomics with optimised shank profile and mono-bloc design
- > Compact glass rod for outstanding illumination of the treatment site
- > Durable instrument construction with stainless-steel external parts
- > Easy change of instruments
- > Thermowasher disinfected and sterilizable
- > 36 months warranty

### SYNEA VISION CONTRA-ANGLE HANDPIECES



1:5, 5-fold spray, head size Ø 9,5mm, press-button chuck system for FG burs Ø 1,6 mm, glass rod

**WK-99 LT**  
[ceramic ball bearings]  
REF 30023000




1:4,5, 5-fold spray, head size Ø 9,1 mm, press-button chuck system for FG burs Ø 1,6 mm, glass rod

**WK-93 LT**  
[ceramic ball bearings]  
REF 30024000



1:1, single spray, head size Ø 9,5 mm, press-button chuck system for contra-angle burs Ø 2,35 mm, glass rod

**WK-56 LT**  
REF 30025000



2:1, single spray, head size Ø 9,5 mm, press-button chuck system for contra-angle burs Ø 2,35 mm, glass rod

**WK-66 LT**  
REF 30026000



10:1, single spray, head size Ø 9,5mm, press-button chuck system for contra-angle burs Ø 2,35 mm, glass rod

**WK-86 LT**  
REF 30027000



### SYNEA VISION STRAIGHT HANDPIECE



1:1, single spray, twist tension chuck system, for handpiece and contra-angle burs Ø 2.35 mm, glass rod

**HK-43 LT**  
REF 30028000




## CONTRA-ANGLE HANDPIECES SHORT EDITION

synea VISION




Short Edition: a significantly lighter contra-angle handpiece system for micromotors with a coupling length up to 23 mm.

- > 7% lighter, 10% shorter
- > Perfect balance and ergonomics of the overall system
- > Innovative, scratch-resistant surface for complete protection
- > Available in the classic range including two different head sizes for the highspeed handpiece
- > Only for micromotors with an ISO coupling length of 23 mm or less
- > Thermowasher disinfected and sterilizable
- > 36 months warranty




1:5, 5-fold spray, head size Ø 9,5mm, press-button chuck system for FG burs Ø 1,6 mm, glass rod

**WK-99 LT S**  
[ceramic ball bearings]  
REF 30215000




1:4,5, 5-fold spray, head size Ø 9,1 mm, press-button chuck system for FG burs Ø 1,6 mm, glass rod

**WK-93 LT S**  
[ceramic ball bearings]  
REF 30216000



1:1, single spray, head size Ø 9,5 mm, press-button chuck system for contra-angle burs Ø 2,35 mm, glass rod

**WK-56 LT S**  
REF 30217000



2:1, single spray, head size Ø 9,5 mm, press-button chuck system for contra-angle burs Ø 2,35 mm, glass rod

**WK-66 LT S**  
REF 30218000



## TURBINE HANDPIECES

synea FUSION



Synea Fusion turbine handpieces – outstanding price-performance ratio. Offering uncompromising user comfort with optimal LED-illumination, Quattro Spray and quiet, vibration-free operation.

- > Largest illumination field with LED+
- > Optimum access, best view and high power with perfect head size to performance ratio
- > Safe cooling with Quattro Spray (4 port spray)
- > Ceramic ball bearings for quiet, vibration-free operation
- > Hygienic head system
- > Thermo washer disinfectable and sterilizable
- > 12 months warranty

### SYNEA FUSION TURBINE HANDPIECES



Quattro Spray, 19 W, head size Ø 11,5 mm, press-button chuck system for FG turbine burs Ø 1,6 mm

**TG-98 L**   
 [Roto Quick connection\*]  
 REF 30001000

**TG-98 L RM**   
 [standard 6-hole connection]  
 REF 30198000

**TG-98**   
 [Roto Quick connection\*]  
 REF 30002000




Quattro Spray, 16 W, head size Ø 10 mm, press-button chuck system for FG turbine burs Ø 1,6 mm

**TG-97 L**   
 [Roto Quick connection\*]  
 REF 30004000

**TG-97 L RM**   
 [standard 6-hole connection]  
 REF 30197000

**TG-97**   
 [Roto Quick connection\*]  
 REF 30005000  
 See Roto Quick couplings

### SYNEA FUSION TURBINE HANDPIECES WITH NSK CONNECTION



Glass rod, Quattro Spray, 19 W, head size Ø 11,5 mm, press-button chuck system for FG turbine burs Ø 1,6 mm

**TG-98 L N**   
 REF 30200000



glass rod, Quattro Spray, 16 W, head size Ø 10 mm, press-button chuck system for FG turbine burs Ø 1,6 mm

**TG-97 L N**   
 REF 30199000

### SYNEA FUSION TURBINE HANDPIECES WITH MULTIFLEX® CONNECTION\*



Glass rod, Four Hole Spray, 19 W, head size Ø 11,5 mm, press-button chuck system for FG turbine burs Ø 1,6 mm






**TG-98 LM**   
 REF 30003000



Glass rod, Quattro Spray, 16 W, head size Ø 10 mm, press-button chuck system for FG turbine burs Ø 1,6 mm

**TG-97 LM**   
 REF 30006000

See W&H coupling with LED light

-  with LED+ light
-  with light
-  without light
-  integrated spray
-  thermo washer disinfectable
-  sterilizable

## CONTRA-ANGLE AND STRAIGHT HANDPIECES

synea FUSION



Synea Fusion – be assured with proven Synea quality at an outstanding price:

- > Long life span and quiet running
- > Optimum access, best view due to small head-sizes
- > Safe cooling
- > Monobloc design for optimum hygiene and ergonomics
- > Compact glass rod for constant light quality
- > Thermowasher disinfected and sterilizable
- > 24 months warranty

### SYNEA FUSION CONTRA-ANGLE HANDPIECES



1:5, Quattro Spray, head size Ø 10 mm, press-button chuck system for FG burs Ø 1,6 mm, glass rod  
Ceramic ball bearings

WG-99 LT  
REF 30007000



WG-99 A  
REF 30008000



2:1, single spray, head size Ø 9,5 mm, press-button chuck system for contra-angle burs Ø 2,35 mm, glass rod

WG-66 LT  
REF 30011000



WG-66 A  
REF 30012000



1:1, Single Spray, head size Ø 9,5 mm, press-button chuck system for contra-angle burs Ø 2,35 mm, glass rod

WG-56 LT  
REF 30009000



WG-56 A  
REF 30010000



### SYNEA FUSION STRAIGHT HANDPIECE




1:1, single spray, twist tension chuck system, for handpiece and contra-angle burs Ø 2.35 mm

HG-43 A  
REF 30013000



## RECIPROCATING CONTRA-ANGLE HANDPIECES

synea FUSION



W&H Profin® contra-angle handpiece for Dentatus-tips:

- > Safe and efficient reworking of excess filler material, cement and bonding
- > Thorough cleaning and polishing of proximal and interdental areas, as well as subgingival surfaces as well as Fine polishing with wooden tips
- > Shaping and adjustment of papillary areas to natural tooth forms
- > Grinding of concave front tooth guiding surfaces
- > Fine processing of grooves and cavities

W&H Intensiv contra-angle handpiece for Intensiv files:

- > Finishing of interproximal and subgingival areas and cavity margins
- > Chamfering of the preparation edges and of cavities
- > Removal of sub gingival plaque and Orthodontic stripping

### SYNEA FUSION RECIPROCATING CONTRA-ANGLE HANDPIECES



Profin® contra-angle handpiece 2:1, single spray, stroke 1.1 mm, for Dentatus tips

WG-67 LT  
REF 30084000



WG-67 A  
REF 30085000



Intensiv contra-angle handpiece 2:1, single spray, stroke 0.9 mm, for Intensiv files

WG-69 LT  
REF 30076000



WG-69 A  
REF 30077000





## TURBINE HANDPIECES

alegra



Discover the Alegra turbines with extremely quiet operation for daily use in your practice. Designed to meet your needs with state-of-the-art technology.

- > Optimal view and illumination of the treatment site thanks to generator technology in the coupling, optimal position of the LED and the small turbine head
- > Extremely quiet operation and long service life due to ceramic ball bearings in turbines
- > Vibration-free, quiet operation thanks to precise manufacturing
- > Perfect spray function thanks to internal water filters
- > Fast removal of tooth mass thanks to extremely high torque
- > Optimal cooling of the treatment site
- > Simple, rapid bur changes
- > Fatigue-free working

### ALEGRA TURBINE HANDPIECES WITH CERAMIC BALL BEARINGS WITHOUT LIGHT



Triple spray, press-button chuck system, max. 16 W head size Ø 12.2 mm, for FG turbine burs Ø 1.6 mm

**TE-98**  
[Roto Quick connection\*]  
REF 30071000

**TE-98 BC**  
[Borden 2-hole connection]  
REF 30072000

**TE-98 RM**  
[Standard 4-hole connection]  
REF 30073000



Triple spray, press-button chuck system, max. 14 W, head size Ø 10.4 mm, for FG turbine burs Ø 1.6 mm

**TE-97**  
[Roto Quick connection\*]  
REF 30068000

**TE-97 BC**  
[Borden 2-hole connection]  
REF 30069000

**TE-97 RM**  
[Standard 4-hole connection]  
REF 30070000

### ALEGRA TURBINE HANDPIECES WITH STEEL BALL BEARINGS WITHOUT LIGHT



Single spray, press-button chuck system, max. 16 W head size Ø 12.2 mm, for FG turbine burs Ø 1.6 mm

**TE-95 BC**  
[Borden 2-hole connection]  
REF 30074000

**TE-95 RM**  
[Standard 4-hole connection]  
REF 30075000

### ALEGRA TURBINE HANDPIECES WITH CERAMIC BALL BEARINGS AND LIGHT



Triple spray, press-button chuck system, head size Ø 12.2 mm, max. 16 W, for FG turbine burs Ø 1.6 mm

**TE-98 LQ**  
REF 30067000



Triple spray, press-button chuck system, head size Ø 10.4 mm, max. 14 W, for FG turbine burs Ø 1.6 mm

**TE-97 LQ**  
REF 30066000

### ROTO QUICK COUPLINGS WITH GENERATOR FOR ALEGRA TURBINE HANDPIECES TE-97/98 LQ



Borden 2-hole connection, generator and water filter with non suck-back stop

**RQ-53**  
REF 10405300



Standard 4-hole connection, generator and water filter with non suck-back stop

**RQ-54**  
REF 10405400

led with LED light

with light

without light

integrated spray

thermo washer disinfectable

sterilizable

See Roto Quick and W&H coupling with LED light

**CONTRA-ANGLE  
AND STRAIGHT  
HANDPIECES**

alegra



Treatment of your patients is very efficient and extremely safe with Alegra. Straight and contra-angle handpieces designed with proven W&H quality.

- > Fatigue-free working thanks to excellent ergonomics
- Thermo washer disinfectable
- > Optimal view of the treatment site due to small head
- > Simple, rapid bur changes with low actuating force thanks to push-button chuck system
- > Optimal cooling of the treatment site thanks to external spray clip (optional)
- > Maximum robustness and durability thanks to high-quality materials

**ALEGRA CONTRA-ANGLE AND STRAIGHT HANDPIECES**



1:1, press-button chuck system, head size Ø 10.1 mm, for contra-angle burs Ø 2.35 mm

**WE-56**  
REF 30116000



1:1, press-button chuck system, head size Ø 10.1 mm, for FG burs Ø 1.6 mm

**WE-57**  
REF 30117000



4:1, press-button chuck system, head size Ø 10.1 mm, for contra-angle burs Ø 2.35 mm

**WE-66**  
REF 30118000



1:1, twist tension chuck system, for handpiece and contra-angle burs Ø 2.35 mm

**HE-43**  
REF 30119000



**ALEGRA CONTRA-ANGLE HANDPIECES WITH GENERATOR**



1:4,5, triple spray, press-button chuck system, head size Ø 10.1 mm, for FG burs Ø 1.6 mm

**WE-99 LED G**  
REF 10259900



4:1, single spray, press-button chuck system, head size Ø 10.1 mm, for contra-angle burs Ø 2.35 mm

**WE-66 LED G**  
REF 10256600



1:1, single spray, press-button chuck system, head size Ø 10.1 mm, for contra-angle burs Ø 2.35 mm

**WE-56 LED G**  
REF 10255600



**SPRAY CLIPS AVAILABLE**

**CONTRA-ANGLE  
AND STRAIGHT  
HANDPIECES**

proxeo



Professional tooth cleaning is an effective method, giving your patients a beaming smile. With Proxeo, W&H offers a compact, easy to handle solution for everything you need, no matter whether it is cleaning, polishing or fluoridating!

**The advantages are clear:**

- > Long life with unique sealing system
- > Best view with small head
- > Efficient treatment with ergonomic design
- > Optimal speeds for atraumatic and efficient polishing/cleaning
- > Easy to clean and hygiene-friendly
- > Thermo washer disinfectable and sterilizable



4:1, contra-angle handpiece system UNIVERSAL with Snap-on adaptor and key for all standard Screw-in cups, brushes and Snap-on adaptor

**WP-64 MU**  
REF 10216443



4:1, contra-angle handpiece for all cups, brushes and Snap-on adaptor with 2,35 mm shank

**WP-66 M**  
REF 30190000



4:1, contra-angle handpiece system YOUNG with key, for YOUNG Screw-in cups and brushes

**WP-64 M**  
REF 10216403



4:1, straight handpiece for disposable contra-angles with plastic shank, system Doriot

**HP-44 M**  
REF 10114403

## ROTO QUICK COUPLINGS

roto quick



### Roto Quick couplings

The Roto Quick system coupling is a multi-faceted system interface, which connects any standardized turbine tubing to any W&H turbine, air motor or air scaler with Roto Quick connection.

- > Fast and safe instrument interchange
- > Free and easy rotation of instruments through 360°
- > Non-retraction valve
- > Sterilizable – for optimum hygiene

### ROTO QUICK COUPLINGS - OPTIC



Standard 6-hole connection

**RQ-24**  
REF 10402400




Standard 6-hole connection,  
adjustable spray

**RQ-34**  
REF 10403400



### W&H COUPLING WITH LED LIGHT FOR KAVO MULTIFLEX® CONNECTION



Standard 6-hole connection,  
adjustable spray

**RM-34 LED**  
REF 10403410



### ROTO QUICK COUPLINGS - NON OPTIC



Borden 2(3)-hole connection

**RQ-03**  
REF 10400300



standard 4-hole connection

**RQ-04**  
REF 10400400



standard 4-hole connection,  
adjustable spray

**RQ-14**  
REF 10401400



## SERVICE OIL F1

f1 service oil



W&H Service Oil F1 is a synthetic oil of maximum purity, which has been specially developed for maintenance of W&H turbine handpieces, straight and contra-angle handpieces, surgical instruments, air scalers and air motors.

- > Longer life for your instruments
- > High material compatibility
- > Sterilization resistant
- > Enormously economical




oil spray can, package of 6 cans


**MD-400 (400 ml, 6 pcs)**  
REF 10940021

Keep spray caps with nozzles  
for oil spray cans in stock  
They can be re-used with  
every new can.


### SPRAY CAPS WITH NOZZLES




for contra-angle and  
straight handpieces  
REF 02038200



for W&H turbine handpieces,  
air motors and air scalers  
with Roto Quick connection  
REF 02229200



for turbine handpieces and air  
motors with fixed connection  
REF 02036100



for products with  
Multiflex® connection  
REF 04396100

## AIR MOTORS



### Air motors AM-25

- > Optimal view and illumination by LED
- > Fatigue-free working thanks to easy 360° rotation
- > Fast removal of tooth structure thanks to high torque

### AIR MOTORS SYSTEM E 5,000 – 25,000 RPM



Standard 4-hole fix connection, for contra-angle and straight handpieces without light

**AM-25 E RM**  
REF 30105000



Borden 2(3)-hole fix connection, for contra-angle and straight handpieces without light

**AM-25 E BC**  
REF 30111000



### AIR MOTORS WITHOUT SPRAY 5,000 – 25,000 RPM



standard 4-hole fix connection, for contra-angle and straight handpieces without light

**AM-25 RM**  
REF 30108000



Borden 2(3)-hole fix connection, for contra-angle and straight handpieces without light

**AM-25 BC**  
REF 30112000



### AIR MOTORS SYSTEM LT / A 5,000 – 25,000 RPM



standard 6-hole connection, for contra-angle and straight handpieces with light

**AM-25 L RM**  
REF 30102000



standard 4-hole fix connection, for contra-angle and straight handpieces without light

**AM-25 A RM**  
REF 30103000



Borden 2(3)-hole fix connection, for contra-angle and straight handpieces without light

**AM-25 A BC**  
REF 30104000



- |                   |   |
|-------------------|---|
| <b>System LT:</b> | with light (light source in motor) and integrated spray |
| <b>System A:</b>  | without light and integrated spray                      |
| <b>System E:</b>  | without light and external spray                        |

- with LED light
- without light
- with light
- integrated spray
- external spray
- without spray
- sterilizable



## PIEZO SCALER

With the Tigon+/Tigon, W&H has developed an exceptional piezo scaler. Its unique features greatly simplify daily work procedures, and patients can relax during treatment as they never could before.

- > Temperature-controlled fluid
- > Perfectly coordinated range of tips
- > Handpiece with 5x Ring LED
- > Large, easy-to-read display
- > Five different treatment programs
- > Three modes for reliable treatment
- > Thermowasher disinfectable and sterilizable: handpiece and tank
- > Multifunctional foot control

### TIGON+ PIEZO SCALER



#### Tigon+, consists of:

- > Control unit (230 V or 115 V models) with 1.8 m supply hose and coupling
- > Handpiece with 5x Ring LED
- > Universal tips set consists of: tip 1x 1U, 1x 3U, 1x 1P, 3x tip changer with torque limitation and tip card
- > Endo/implant tip set consists of: tip 1x 1I, 1x 1E, 1x tip changer with torque limitation, 1x open-ended spanner and tip card
- > Foot control S-N1
- > 400 ml tank
- > Coolant hose with coolant filter
- > Tray
- > Handpiece support
- > Nozzle cleaner

**PA-123**  
REF 13931000

### TIGON PIEZO SCALER




#### Tigon, consists of:

- > Control unit (230 V or 115 V models) with 1.8 m supply hose and coupling
- > Handpiece with 5x Ring LED
- > Universal tips set consists of: tip 1x 1U, 1x 3U, 1x 1P, 3x tip changer with torque limitation and tip card
- > Foot control C-NF
- > 400 ml tank
- > Tray
- > Handpiece support
- > Nozzle cleaner


**PB-323**  
REF 13923000

### PIEZO TIPS SETS FOR TIGON+/TIGON



- Consists of:
- > tip 1U, 3U, 1P (1 pc)
  - > tip changer with torque limitation (3 pcs)
  - > tip card

**Universal tip set**  
REF 05851800\*



- Consists of:
- > tip 2Pr, 2PI, 1P (1 pc)
  - > tip changer with torque limitation (3 pcs)
  - > tip card

**Perio tip set**  
REF 06370200



- Consists of:
- > tip 1I (3 pcs)
  - > tip changer with torque limitation (3 pcs)

**Implant clean tip set**  
REF 06370500



### PIEZO TIPS SETS FOR TIGON+/TIGON



- Consists of:
- > tip 1E, 2E, 1R (1pc)
  - > spanner (1 pc) for tip 1E, 2E
  - > tip changer with torque limitation (1 pc) for tip 1R
  - > tip card

**Endo tip set**  
REF 06370400



- Consists of:
- > tip 1E (3 pcs)
  - > spanner (1 pc)
  - > tip card

**Endo tip set**  
REF 06370600



- Consists of:
- > tip 1E, 1I (1 pc)
  - > spanner (1 pc) for tip 1E
  - > tip changer with torque limitation (1 pc) for tip 1I
  - > tip card

**Endo / implant tip set**  
REF 06382800





## AIR SCALER

proxeo

Fast, gentle, powerful. The air-driven scaler with the integrated spray and two large optic outlets is ideal for the removal of plaque and calculus as well as cleaning root surfaces. The result: your patients feel better in general and more self-confident thanks to their clean teeth surfaces.

- > Stepless regulation
- > Large optic outlets
- > Large tip range
- > Integrated spray
- > First-class hygiene
- > Thermo washer desinfectable and sterilizable

### PROXEO AIR SCALER WITH ROTO QUICK CONNECTION



**ZA-55 L**  
REF 10035500



with power regulation and tips 1AU, 2AU, 3AU and 1AP; tip changer, nozzle cleaner, small cleaning brush, spare bulb, allen key and tip card



**ZA-55**  
REF 10035501



with tips 1AU, 2AU, 3AU and 1AP; tip changer, nozzle cleaner, small cleaning brush and tip card

### PROXEO AIR SCALER WITH FIXED CONNECTION



**ZE-55 BC**  
REF 11015500



with tip 1AU (2 pcs), open-ended spanner, nozzle cleaner, small cleaning brush and tip card



**ZE-55 RM**  
REF 11015501



with tip 1AU (2 pcs), open-ended spanner, nozzle cleaner, small cleaning brush and tip card

### PROXEO AIR SCALER WITH SIRONA® CONNECTION



**ZA-55 LS**  
REF 30227000



with power regulation and tips 1AU, 2AU, 3AU and 1AP; tip changer, nozzle cleaner, small cleaning brush and tip card

### PROXEO AIR SCALER WITH MULTIFLEX® CONNECTION\*



**ZA-55 LM**  
REF 10035502



with power regulation and tips 1AU, 2AU, 3AU and 1AP; tip changer, nozzle cleaner, small cleaning brush and tip card



**ZA-55 M**  
REF 10035503



with tips 1AU, 2AU, 3AU and 1AP; tip changer, nozzle cleaner, small cleaning brush and tip card

### TIPS FOR PROXEO AIR SCALER



**Tip 1AU with tip card**  
REF 07159700



Universally applicable, also for curved molar surfaces with lingual inclination



**Tip 2AU with tip card**  
REF 07159800



For anterior teeth of the lower jaw, lingual surfaces and in cases of stubborn adhesion of thick, hard collars of dental calculus



**Tip 3AU with tip card**  
REF 07159900



For removing plaque from the neck of the tooth and removing nicotine stains



**Tip 1AP with tip card**  
REF 07009100



Removal of subgingival deposits, especially suitable for deep periodontal pockets



**Tip 2APr with tip card**  
REF 05254400



Tip for removal of subgingival deposits. Right-curved and leftcurved models available



**Tip 2API with tip card**  
REF 05153300



Tip for removal of subgingival deposits. Right-curved and left curved models available



**Tip 3AP with tip card**  
REF 05280200



With a diamond coating on the instrument tip for periodontal debridement of furcations and concavities.

W&H adaptor kit for cleaning tips in thermo washer disinfectors available